

Modular Coil Geometry for Case 0711a1 and 0806a1


Case 0806a1:

Min plasma-coil dist = 19-cm

Min coil-coil dist = 11.3-cm

Min NB-coil dist = 25.2-cm


Min rad of curvature = 9.1-cm



Segmentation Options for 18 Modular Coils (0806a1)


- PVR segmentation (Shell and TF) is 8.6, 25.7, 42.9, 60-deg

18 TF Coils



Shell segments same as TF

24 TF Coils



Shell segments at 0,21,36,60-deg